

SILICONES - CERTIFICATE OF ANALYSIS

B220005 ROYAL BITE 2x50ML

MANUFACTURED DATE: 2023-10-04	BATCH NUMBER: 0000431734	EXPIRY DATE: 2025-10-02
--------------------------------------	---------------------------------	--------------------------------

Test Conditions	Test Results	Specifications	Test method
Room temperature (°C)	22,4	[22,0-24,0]	2100-WI QC-001
Relative Humidity (%RH)	52,1	[40,0-60,0]	2100-WI QC-001
Base temperature (°C)	22,9	[22,0-24,0]	2100-AM QC-023
Catalyst temperature (°C)	22,7	[22,0-24,0]	2100-AM QC-023
Mixing time (s)	Compliant	0 (Automatic mix	

Features	Test Results	Specifications at release ¹	Warranty specifications ²	Test method
Colour base	Compliant	Coral	Coral	2100-AM QC-004
Colour catalyst	Compliant	White	White	2100-AM QC-004
Visual check	Compliant	Homogeneous*	Homogeneous*	2100-AM QC-004
Mixed material	Compliant	Homogeneous*	Homogeneous*	2100-AM QC-004
Hardness 8' (23°C) (ShA)	84,0	[83,0-93,0]	83 - 93	2100-AM QC-013
Total working time 35°C (s)	34	[30-50]	30 - 50	2100-AM QC-023
Setting time 35°C (s)	230	[160-260]	160 - 300	2100-AM QC-023

Visual Inspection Report

Finished product inspection result: COMPLIANT.

Badia Polesine, 2023-10-17 QUALITY CONTROL LABORATORY

1. Values are referred to the date when the product is released.
2. Values are referred to the features of the product up to the expiry date.

*Homogeneous and free of extraneous matter

This is an electronic document, the signature is not necessary.