

## SILICONES - CERTIFICATE OF ANALYSIS

**B220001 ROYAL DENT PUTTY SOFT FAST SET 300+300ML**

<b>MANUFACTURED DATE: 2023-10-03</b>	<b>BATCH NUMBER: 0000431712</b>	<b>EXPIRY DATE: 2025-10-01</b>
--------------------------------------	---------------------------------	--------------------------------

Test Conditions	Test Results	Specifications	Test method
Room temperature (°C)	22,4	[22,0-24,0]	2100-WI QC-001
Relative Humidity (%RH)	51,2	[40,0-60,0]	2100-WI QC-001
Base temperature (°C)	22,6	[22,0-24,0]	2100-AM QC-023
Catalyst temperature (°C)	22,7	[22,0-24,0]	2100-AM QC-023
Mixing time (s)	Compliant	30	

Features	Test Results	Specifications at release <sup>1</sup>	Warranty specifications <sup>2</sup>	Test method
Colour base	Compliant	Turquoise	Turquoise	2100-AM QC-004
Colour catalyst	Compliant	White	White	2100-AM QC-004
Visual check	Compliant	Homogeneous*	Homogeneous*	2100-AM QC-004
Mixed material	Compliant	Homogeneous*	Homogeneous*	2100-AM QC-004
Hardness 8' (35°C) (ShA)	61	[57-63]	58 - 65	2100-AM QC-013
Total working time 35°C (s)	77	[70-100]	65 - 100	2100-AM QC-023
Setting time 35°C (s)	159	[135-170]	120 - 160	2100-AM QC-023

### Visual Inspection Report

Finished product inspection result: COMPLIANT.

Badia Polesine, 2023-10-17      QUALITY CONTROL LABORATORY

1. Values are referred to the date when the product is released.
2. Values are referred to the features of the product up to the expiry date.

\*Homogeneous and free of extraneous matter

This is an electronic document, the signature is not necessary.